

SUPER-BOND

UNIVERSAL CERAMIC PRIMER


Designed for All Types of Ceramic and Zirconia
Assured Bonding Performance



A priming agent for ceramic surface. Application of Super-Bond Universal Ceramic Primer on a ceramic surface before bonding with a dental adhesive Resin cement will enhance the bond strength and durability.

Bond Strength with Zirconia

Tensile Bond Strength with Different Zirconia* (Mpa)




* Tested after thermocycled 10,000 times (5-55°C).
Resin Cement: Super-Bond

Assures to bond zirconia firmly to tooth using Super-Bond

Bond Strength with Porcelain

Tensile Bond Strength with Different Porcelain* (Mpa)



* Tested after thermocycled 10,000 times (5-55°C).
Resin Cement: Super-Bond

Guarantees a reliable adhesion to porcelain.

Clinical Case

A Zirconia Bridge



Application of Super-Bond Universal Ceramic Primer



Dispense Liquid A and Liquid B, mix and apply to the inner surface of the zirconia bridge

Application of Resin Cement



Use Super-Bond to bond the bridge

After setting



Contents of the Kit

1. Liquid A - 5ml
2. Liquid B - 5ml