

HYBRID ROOT SEAL

SELF ETCHING RESIN-BASED DUAL-CURED ROOT CANAL SEALER

"Don't just fill, seal your root canals

Hybrid Root SEAL is a Self-etching, Resin Based, dual cured root canal sealer which hybridizes both canal walls & Gutta Percha to prevent leakage.

Why a Dentist should buy Hybrid Root SEAL?

- ▣ Extremely easy to use.
- ▣ Dual Cure.
- ▣ Forms Hybrid Layer with Canal Wall & Gutta Percha.
- ▣ Hydrophilic Polymerization Initiator.
- ▣ Long lasting results.
- ▣ Freedom to use it in all canal types.
- ▣ Compatible with all cold obturation techniques.
- ▣ Radiopaque.

Get unexpected results with Hybrid Root SEAL

- ▣ Highly Bio-compatible.
- ▣ Prevents leakage during provisional period.
- ▣ Excellent sealing with only one cone.
- ▣ Polymerization starts even on Humid surface of root canal.



Where does company recommend Hybrid Root SEAL

Root Canal Sealer for all types of Root Canals compatible with any point.

Technical data

1. Micro Tensile Bond Strength

A. EDTA/NaClO (60 sec) - 14.5 + 2.8 MPa

B. None - 17.7 + 4.2 MPa


2. Sealability

A. Immediately - 95-100%

B. Durability - 95-100%

Instructions for use

How to use 'Hybrid Root SEAL'



Packaging

1. Powder 5 gm
2. Liquid 4.5 ml
3. Spatula
4. Measuring Spoon
5. Mixing Pad

