

# METAFIL CX

## RESIN-BASED DENTAL RESTORATIVE MATERIAL

Excellent Mechanical Property  
Stress-free Non-sticky Handling



## Outstanding Polishability



Metafil CX is a microfill, natural-looking resin composite restorative using a unique pre polymerized organic filler TMT, which is very reactive with the matrix. This feature gives the finished restoration greater resilience and wear resistance. Metafil CX has excellent adaption and non-sticky handling. It is easy to place beautiful restorations even in the most difficult scenarios. Shades are precisely matched to the VITA shade guide making shade matching both simple and accurate. Metafil CX is super-polishable, it gives outstanding esthetic after quick polishing.

## Where does company recommend Metafil CX

"Make your Anterior Restorations Invisible"

Primarily Indicated for ANTERIOR & CERVICAL restorations:

- ▣ Direct Veneers.
- ▣ Diastema Closure.
- ▣ Class III, Class IV & Class V cavity fillings.
- ▣ Class I & Class II- it's radiolucent !

## Technical data

Compressive Strength	-	420 MPa
Flexural Modulus	-	3.5 GPa
Refractive Index	-	1.48

## Typical Applications

Fractured incisors



Shade matching on teeth



Restored with Metafil CX



After polishing



Superior polishability of Metafil CX will guarantee highly esthetic anterior restoration

Pre-operative labial view



Pre-operative incisal view



Labial view after restoration



Incisal view after restoration



Metal CX is very retentive to enamel, and it is suitable for closing diastemas.

Pre-operative labial view



After prep and application of bonding agent



Restoration with Metafil CX and polish



After completion



After quick polish, Metafil CX is highly esthetic and ideal for anterior restoration

## Shades



## Packaging

Content of each syringe-3.5 gm